

Liquid Urea Nitrogen (BUN) Reagent Set

Intended Use

For the quantitative determination of urea nitrogen in serum using the Mindray BS-200 analyzer. For *in vitro* diagnostic use only.



Clinical Significance

Determination of urea nitrogen in serum is widely used as a screening test for renal function. When used in conjunction with the determination of creatinine in serum it is helpful in the differential diagnosis of the three types of azotemia; pre-renal, renal and post-renal.¹

Method History

Urea has been determined by the direct method² where urea condenses with diacetyl to form a chromagen and an indirect method where ammonia is measured as a product of urease action on urea.³ The liberated ammonia has been measured using Nessler's reagent⁴ and by the Berthelot reaction.⁵ Talke and Schubert introduced a totally enzymatic procedure in 1965 utilizing urease and glutamate dehydrogenase.⁶ The present procedure is based on a modification of their method.

Principle



Urea is hydrolyzed by urease to produce ammonia and carbon dioxide. The liberated ammonia reacts with α -ketoglutarate in the presence of NADH to yield glutamate. An equimolar quantity of NADH undergoes oxidation during the reaction resulting in a decrease in absorbance that is directly proportional to the urea nitrogen concentration in the sample.

Reagent Composition

Working reagent concentrations: Urease (Jack Bean) >15,000 U/L, GLDH (Bovine) >200 U/L, ADP >0.6 mM, α -Ketoglutarate 3.4 mM, NADH >0.28 mM, Buffer, stabilizers, Sodium Azide (0.28%) as preservative.

Reagent Preparation

The reagents are ready to use.

Reagent Storage

Store R1 and R2 reagents at 2-8°C. The reagents are stable until the expiration date appearing on the label when stored as directed.

Reagent Deterioration

The reagent should not be used if the working reagent has a reagent blank absorbance less than 1.0 at 340 nm.

Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Avoid ingestion of reagent as toxicity has not yet been determined.
3. Reagents contain sodium azide (0.28%) as preservative. Sodium azide may react with copper or lead plumbing to form explosive metal azides. Upon disposal flush with large amounts of water.

4. All specimens should be handled in accordance with good laboratory practices using appropriate precautions as described in the CDC/NIH Manual "Biosafety in Microbiological and Biomedical Laboratories," 2nd ed., 1988, HHS Publication No. (CDC) 88-8395.

Specimen Collection and Storage

1. Serum is recommended.
2. Plasma containing anticoagulants should not be used.
3. All material coming in contact with the sample must be free of ammonia and heavy metals.⁷
4. Urea in serum is reported stable for seventy-two hours refrigerated at 2-8°C. Unrefrigerated sera should be used within eight hours.
5. Specimen collection should be carried out in accordance with NCCLS M29-T2.⁸ No method can offer complete assurance that human blood samples will not transmit infection. Therefore, all blood samples should be considered potentially infectious.

Interferences

1. Urease action is inhibited by fluoride.
2. Samples with abnormal ammonia levels give falsely elevated BUN results.
3. Bilirubin to the level of 20 mg/dl was found to exhibit negligible interference (<2%) in this assay.
4. Hemoglobin to the level of 200 mg/dl was found to exhibit negligible interference (<5%) in this assay.
NOTE: The BUN level was 46.0 mg/dl for the Bilirubin study and 46.3 mg/dl for the Hemoglobin study.
5. For a comprehensive review of drug interference see Young, et al.⁹

Materials Provided

Urea Nitrogen Enzyme Reagent (R1)
Urea Nitrogen Coenzyme Reagent (R2)

Materials Required but not Provided

1. Mindray BS-200 Analyzer
2. BS-200 Operation manual
3. Chemistry Calibrator, catalog number C7506-50
4. Chemistry control, catalog number C7592-100

BS-200 Test Parameters

Test:	BUN	R1:	300
No.:	008	R2:	75
Full Name:	Urea Nitrogen	Sample Volume:	3
Standard No.:		R1 Blank:	
Reac. Type:	Fixed-time	Mixed Rgt. Blank:	
Pri. Wave:	340nm	Linearity Range:	0 - 150
Sec. Wave.:	670nm	Linearity Limit:	
Direction:	Decrease	Substrate Limit:	
Reac. Time:	2 / 7	Factor:	
		Compensate: Slope 1.0 Intercept 0	
Incuba. Time:	3	<input type="checkbox"/> Prozone check	
Unit:	mg/dl	q1: q2: q3: q4:	
Precision:	Integer	PC: Abs:	

Liquid Urea Nitrogen (BUN) Reagent Set

Calibration Parameters

Rule:	Two-point linear	Calibrator 1:	Deionized Water
Sensitivity:		Calibrator 2:	Chem Cal
Replicates:	2	Calibrator 3:	
Interval (day):		Calibrator 4:	
Difference Limit:		Calibrator 5:	
SD:		Calibrator 6:	
Blank Response:			
Error Limit:			
Coefficient:	0		

Limitations

Samples with values above 150 mg/dl should be diluted with 0.9% saline 1:1, re-assayed and the results multiplied by two.

Calibration

Use an NIST-traceable serum calibrator. The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the procedure should be recalibrated.

Calculation (Example)

$(A_1 - A_2)$ = Absorbance change between readings

$\frac{(A_1 - A_2) \text{ unknown}}{(A_1 - A_2) \text{ standard}} \times \text{concentration of standard} = \text{BUN (mg/dl)}$

Example: If the unknown had an $A_1 = 1.5$ and $A_2 = 1.0$,
the standard $A_1 = 1.5$ and $A_2 = 0.9$ and
the concentration of the standard = 20 mg/dl then:

$$\frac{(1.5 - 1.0)}{(1.5 - 0.9)} = \frac{0.5}{0.6} \times 20 = 17 \text{ mg/dl}$$

NOTE: To obtain results in SI units multiply by 10 to convert dl to liters and divide by 28, the molecular weight of nitrogen.

Example: $17 \text{ mg/dl} \times 10/28 = 6.06 \text{ mmol/L}$.

To convert mg/dl Urea Nitrogen to mmol Urea/L, multiply the mg/dl Urea Nitrogen value by 0.357.

To convert mg/dl Urea Nitrogen to mg/dl Urea, multiply the mg/dl Urea Nitrogen value by 2.14.

Quality Control

The validity of the reaction should be monitored by use of the control sera with known normal and abnormal BUN values. These controls should be run at least with every working shift in which urea nitrogen assays are performed. It is recommended that each laboratory establish its own frequency of control determination. Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

Expected Values

7-18 mg/dl⁷

It is strongly recommended that each laboratory establish its own reference range.

Performance

- Assay Range: 0-150 mg/dl. Samples that exceed 150 mg/dl should be diluted with an equal volume of saline and re-assayed. Multiply the result by two.
- Comparison: A study was performed between the Mindray BS-200 and a similar analyzer using this method, resulting in a correlation coefficient of 0.986 and the regression equation of $y = 0.95x + 0.6$.
- Precision: Precision studies were performed using the Mindray BS-200 analyzer following a modification of the guidelines which are contained in NCCLS document EP5-T2.¹⁰

Within Run			Day to Day		
Mean	S.D.	C.V.%	Mean	S.D.	C.V.%
15.6	0.5	3.2	14.1	0.8	5.7
55.3	1.2	2.1	51.1	1.8	3.5


- Sensitivity: The sensitivity for the Liquid BUN reagent was investigated by reading the change in absorbance at 340 nm for a saline sample, and serum samples with known concentrations. Ten replicates of each sample were performed. The results of this investigation indicated that, on the analyzer used, the Liquid BUN reagent showed little or no drift on a zero sample. Under the reaction conditions described, 1mg/dl of BUN gives an absorbance of 0.003.


References


- Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia W.B. Saunders (1976).
- Fearon, W.R., Biochem J. 331:902 (1939).
- Marshall, E.K., Jr., J. Biol. Chem. 15:487 (1913).
- Gentzkow, C.J., J. Biol. Chem. 143:531 (1952).
- Fawcett, J.K., Scott, J.E., J. Clin. Path. 13:156 (1960).
- Talke, H., Schubert, G.E., Klin. Wschr. 43:174 (1965).
- Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia W.B. Saunders, p991 (1976).
- NCCLS document "Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue", 2nd Ed. (1991).
- Young, D.S., et al, Clin. Chem. 21:1D (1975).
- NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Ed. (1992).

Rev. 10/15

M803-BUN600-01

 Use by (YYYY-MM)

 Temperature limitation

 Lot and batch code


 Consult instructions for use

 Catalog number

 CE mark

 Manufacturer

 Authorized representative in the European Community

 In vitro diagnostic medical device