

autoLDL[™] Cholesterol Reagent Set

Intended Use

For the direct quantitative determination of low density lipoprotein cholesterol (LDL-C) in human serum or plasma using the Mindray BS-200 analyzer. For *in vitro* diagnostic use only.

Summary

Plasma lipoproteins are spherical particles that contain varying amounts of cholesterol, triglycerides, phospholipids, and proteins. The phospholipid, free cholesterol and protein constitute the outer surface of the lipoprotein particle, the inner core contains mostly esterified cholesterol and triglycerides. These particles serve to solubilize and transport cholesterol and triglycerides in the bloodstream.

The relative proportions of protein and lipid determine the density of these plasma lipoproteins and provide a basis for their classification.¹ The classes are: very low density lipoproteins (VLDL), low density lipoproteins (LDL), and high density lipoprotein (HDL). Numerous clinical studies have shown that the different lipoprotein classes have varied effects.²⁻⁴ The studies all point to LDL cholesterol as the key factor in the pathogenesis of artherosclerosis and coronary artery disease (CAD),²⁻⁸ while HDL cholesterol has often been observed to have a protective effect. Even within the normal range of total cholesterol concentrations, an increase in LDL cholesterol can occur with an associated risk for CAD.⁴

Over the years a variety of methods have been employed for the determination, or estimation, of LDL cholesterol. The Friedewald equation, in a variety of forms, has been most frequently used for the estimation of LDL cholesterol. However, its usefulness is limited and its accuracy has been questioned. Determination of LDL cholesterol by beta-quantification is recognized as the reference method, but the procedure is so cumbersome relatively few laboratories use this method. A recent method using immunoseparation has become popular. However, this method is still requires sample pre-treatment prior to cholesterol determination, making it unsuitable for full automation of LDL cholesterol in a two part, liquid stable reagent that is easily adapted to most automated chemistry analyzers.

Reagent Composition

Components	Appearance	Ingredients
Reagent 1	Liquid	MES Buffer (pH 6.3)
		Detergent 1
		Cholesterol esterase
		Cholesterol oxidase
		Peroxidase
		4-aminoantipyrine
		Ascorbic acid oxidase
		Preservative
Components	Appearance	Ingredients
Reagent 2	Liquid	MES Buffer (pH 6.3)
		Detergent 2
		N,N-bis (4-sulfhobutyl)-
		m-Toluidine-disodium
		(DSBmT)
		Preservative

Cholesterol Oxidase from Nocardia sp., Cholesterol Esterase from Pseudomonas sp., Peroxidase from Horseradish, Ascorbic Acid Oxidase from Cucurbita sp.

Principle

The autoLDL[™] Cholesterol Reagent is a two-part, liquid stable method for directly measuring LDL-C levels in serum or plasma. The method depends

on the properties of a unique detergent which eliminates the need for any off-line pre-treatment or centrifugation steps. This detergent (Reagent 1) solubilizes only the non-LDL lipoprotein particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non-color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL-C in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

HDL, VLDL, Chylomicrons	> Solubilized HDL, VLDL, Chylomicrons	>Cons (umed HDL, VLDL, Chylomicrons (No color)
Detergent/Rgt1 LDL	≻ Non-solubilized LDL Cholesterol	Detergent/Rgt2	≻ Solubilized LDL Cholesterol
Solubilized LDL Cholesterol	Cholesterol Esterase Cholesterol Oxidase	> + H	202
H ₂ O ₂ + DSBmT + 4-AA	Peroxidase	≻ Color Dev (Measured Bichro at 546 &	relopment omatically 660nm)

Reagent Preparation

Reagent 1: Reagent 1 is ready to use. Reagent 2: Reagent 2 is ready to use.

Reagent Storage and Stability

All reagents are stable until the expiration date on the label when stored at 2 to 8° C.

Precautions

- 1. Reagent is intended for *in vitro* diagnostic use only.
- 2. Do not pipette by mouth.
- All specimens used in this test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.
- 4. Do not use the reagents beyond the expiration date printed on the kit label.

Specimen Collection and Storage

Serum, EDTA-treated or heparinized plasma are the recommended specimens. Patients are not required to fast prior to blood collection.

- Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).¹⁰
- Plasma: Specimens may be collected in EDTA or heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).¹⁰

If not analyzed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens must be stored for more than 5 days, they may be frozen at -80° C.

Interferences

All interference studies were conducted according to the procedures recommended in NCCLS guideline No. EP7-P for interference testing in clinical

chemistry.¹² Hemoglobin at levels up to 400 mg/dl, Bilirubin at levels up to 20 mg/dl and Triglycerides to 1500 mg/dl were found to exhibit negligible interference (<5%) on this method. Samples with levels of interfering substances higher than the upper limits should be diluted with physiological saline before assaying. Multiply the result obtained from the manual dilution by the appropriate dilution factor. For a comprehensive review of drug interference on serum LDL cholesterol levels see Young et al.¹³

Materials Provided

autoLDL [™] Cholesterol Reagent Set		
Reagent 1	3 x 40mL	
Reagent 2	3 x 14mL	

Materials Required but not Provided

- 1. Pointe Scientific autoHDL/LDL[™] Calibrator, Cat. No. H7545-CAL
- 2. Mindray BS-200 Analyzer
- 3. BS-200 Operation manual
- 4. Pointe Lipid controls, catalog number L7580-18

Procedure

Below is general example of the autoLDL[™] test procedure for an automated analyzer. All analyzer applications should be validated in accordance with NCEP and CLIA recommendations.¹⁰ For assistance with applications on automated analyzers, please contact Pointe Scientific's Technical Service Department at (800)445-9853.

Sample + 3ul	Reagent 1 300ul	37°C ≻ 5min.	Reagent 2 100ul	37°C ≻ 5min.	Measurement (Absorb. Difference between 660nm & 546nm)
					I
					I
					M

LDL-C Result

Mindray BS-200 Test Parameters

Willing DS-200	<u> Test Faraineter</u>	3	
Test:	LDL	R1:	180
No.:	024	R2:	60
Full Name:	autoLDL	Sample Volume:	3
Standard No.:		R1 Blank:	
Reaction Type:	End-point	Mixed Rgt. Blank:	
Pri. Wave:	546nm	Linearity Range:	0 - 370
Sec. Wave:	670nm	Linearity Limit:	
Direction:	Increase	Substrate Limit:	
Reac. Time:	0 / 18	Factor:	
		Compensate: Slope 1.0	Intercept: 0
Incuba. Time:	18	Prozone check	
Unit:	mg/dl	q1: q2: q3: q4:	
Precision [.]	Integer	PC: Abs	

Calibration Parameters

Rule:	Two-point linear	Calibrator 1:	Deionized Water	
Sensitivity:		Calibrator 2:	HDL/LDL Cal	
Replicates:	2	Calibrator 3:		
Interval (day	y):	Calibrator 4:		
Difference L	.imit:	Calibrator 5:		
SD:		Calibrator 6:		
Blank Resp	onse:			
Error Limit:				
Coefficient:	0			

Limitations

- 1. Anticoagulants containing citrate should not be used.
- 2. Protect the reagents from direct sunlight.
- 3. Samples with values greater than 370 mg/dl on the Mindray BS-200 must be diluted 1:1 with saline and re-assayed. Multiply the result by two.

Calibration

The autoHDL/LDL[™] Cholesterol Calibrator is required for calibration. The values of the calibrator were assigned by procedures traceable to the National Reference System for Cholesterol (NRS/CHOL). Refer to autoHDL/LDL[™] Cholesterol Calibrator package insert for instructions. If control results are found to be out of range, the procedure should be re-calibrated.

Quality Control

Reliability of test results should be routinely monitored with control materials that reasonably emulate the performance of patient specimens.¹⁰ Quality control materials are intended for use only as monitors of accuracy and precision. The recovery of control values within the appropriate range should be the criteria used in evaluation of future assay performance. Controls should be run with every working shift in which LDL-C assays are performed. It is recommended that each laboratory establish its own frequency of control determination. Quality control requirements should be determined in conformance with local, state, and/or Federal regulations or accreditation requirements.

Results

To convert from conventional units to S.I. units, multiply the conventional units by 0.02586.

Example: mg/dL x 0.02586 = mmol/L LDL-C

Expected Values

The following NCEP recommendations for patient classifications are suggested for the prevention and management of coronary heart disease.⁸

LDL Cholesterol	<u>Classifications</u>
<130mg/dl (3.36mmol/L) 130-159m/dl (3.36-4.11mmol/L)	Borderline High Risk
160mg/dl (4.14mmol/L)	High Risk

It is highly recommended that each laboratory establish its own range of expected values.

Specific Performance Characteristics Data Generated on BS200

Assay Range: 0-370 mg/dl

Accuracy: Studies comparing the Liquid autoLDL[™] Cholesterol Reagent method used on the Mindray BS200 and a similar analyzer yielded the following results:

Method	autoLDL™	
	Cholesterol	
Ν	30	
Mean LDL Cholesterol	106	
Range (mg/dl)	50-159	
Standard Deviation (mg/dl)	28	
Regression Analysis	Y=1.017x + 1.7	
Correlation Coefficient	R=0.990	



autoLDL[™] Cholesterol Reagent Set

Precision:

Within-Day precision for the autoLDL[™] Cholesterol Reagent was determined following a modification of NCCLS document EP5-T2¹⁷ using the Mindray BS200. Within-Day precision studies produced the following results:

Sample	LOW	HIGH
Ν	20	20
Mean LDL Cholesterol (mg/dl)	28	146
Standard Deviation (mg/dl)	0.5	2.5
Coefficient of Variation (%)	1.7	1.7

Day-to-Day precision was also determined following a modification of NCCLS document EP5-T2.¹⁷ Day-to-Day precision studies run on the Mindray BS200 produced the following results:

Sample	LOW	HIGH
Ν	20	20
Mean LDL Cholesterol (mg/dl)	26	147
Standard Deviation (mg/dl)	0.9	2.9
Coefficient of Variation (%)	3.3	2.0

Sensitivity: 2SD limit of detection (95% Conf) = 0.627 mg/dl.

References

- 1. Gotto, A.M., Lipoprotein Metabolism and the etiology of Hyperlipidemia, Hospital practice, 23:Suppl. 1,4 (1988).
- Crouse, J.R., et al., Studies of Low Density Lipoprotein Molecular Weight in Human Beings with Coronary Artery Disease, J. Lipid Res., 26:566 (1985).
- Badimon, J.J., Badimon L., Fuester V., Regression of Athroscierotic Lesions by High-density lipoprotein Plasma fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85:1234-41 (1990).
- Castelli, W.P., et al., Cholesterol and other Lipids in coronary heart disease, Circulation, 55-767 (1977).
- 5. Barr, D.P., Russ, E.M, Elder, H.A., Protein-Lipid Relationships in Human Plasma, Am. J. Med. 11:480 (1951).
- Gordon, T., et al, High Density Lipoprotein as a Protective Factor Against Coronary Heart Disease, Am. J. Med., 62:707 (1977).
- 7. William, P., Robinson, D., Baily A., High Density Lipoprotein and Coronary Risk Factor, Lancet, 1:72 (1979).
- Kannel, W.B., Castelli W.P., Gordon, T., Cholesterol in the Prediction of Artheriosclerotic Disease; New Perspectives Based on the Framingham Study, Am. Intern. Med., 90:85 (1979).
- 9. National Institutes on Health Publication no. 93-3095, September 1993.
- Warnick, G. Russell, Wood Peter D., National Cholesterol Education Program Recommendations for Measurement of High Density Lipoprotein Cholesterol; Executive Summary, Clinical Chemistry, Vol. 41, No. 10, 1995.
- Grundy, S.M., et al, Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II) JAMA 1993, 269:23,3015-3023.
- National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 4, No. 8, June 1984.
- Young, D.S. Effects of Drugs on clinical Laboratory Tests, 3rd ed., AACC Press, Washington, D.C., 1990, 3-104 thru 3-106.

- 14. Tietz, N.W., Clinical Guide to Laboratory Tests, W.B. Saunders Co., Philadelphia, 1986, p. 256.
- Carey, R., Gerber, C.C., Evaluation of Methods. In Kaplan LA, Pesce, A.J., eds. Clinical Chemistry: theory, analysis and correlation. Third Edition. St. Louis: The CV Mosby Company.
- Westgard, J.O., Carey, R.N., Wold, S., Criteria for judging precision and accuracy in method development and evaluation. Clinical Chemistry 1974:20:825-833.
- NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices" 2nd Ed. 1992.

Symbol Key

Use by (YYYY-MM)	LOT Lot and batch code	
REF Catalog number	Manufacturer	
IVD In vitro diagnostic medical device	X Temperature limitation	
Consult instructions for use	CE mark	
EC REP Authorized representative in the European Community		

Rev: 12/14 M803-LDL600-01

5449 Research Dr. • Canton MI 48188 • USA • Phone: 734-487-8300 • Toll Free: 800-445-9853 • Fax: 734-483-1592 • www.medtestdx.com